

HI-Quik Dengue NS1 Rapid Card Test

For *In-Vitro* Diagnostic Use Only

Store at 4°C to 30°C

INTENDED USE

Dengue NS1 test is an immunochromatographic assay for the qualitative Detection of NS1 Antigen of Dengue virus in human serum/plasma.

PRINCIPLE

Dengue NS1 Antigen Test: After **addition of the serum or plasma sample** to the sample well of the device containing a test strip, the sample moves on to the gold conjugate pad containing colloidal gold particles conjugated with Dengue NS1 antigen specific antibodies and rabbit IgG. If the sample contains detectable levels of the Dengue NS1 antigens it reacts with the gold conjugated Dengue NS1 antibodies to form a complex. This complex along with unbound gold particles move on nitrocellulose membrane. The complex reacts with Dengue NS1 antibodies coated on nitrocellulose membrane at test side to form a colored band (Test Line). The unbound complex, unbound gold conjugate particles and the rabbit IgG conjugated colloidal gold particles move further to the goat anti-rabbit IgG coated control area to form a colored band (Control line). The appearance of test line and control line in respective area indicates the positive result. Appearance of only control line indicates a negative result. The control line acts as a procedural control. Control line should always appear if the test is performed as per the procedure and reagents are working properly.

Please read all the information in this package insert before performing the test.

- ✓ Do not use after the expiration date printed on the foil pouch.
- ✓ Store in the sealed pouch in a dry place at 4°C to 30°C. Do not freeze.
- ✓ Do not use if pouch is torn or damaged.
- ✓ Do not open the foil pouch until you are ready to start the test.
- ✓ Do not reuse the test device.
- ✓ Dispose of hygienically in domestic waste.
- ✓ Do not touch the membrane.
- ✓ Keep out of the reach of children.
- ✓ For in vitro diagnostic use. Not to be taken internally.
- ✓ Do not eat the desiccant in the package.

CONTENTS OF KIT

1. Test Device & Silica Gel in pouch
2. Package Insert
3. Plastic dropper for Dengue NS1

MATERIALS NEEDED BUT NOT PROVIDED

1. Timer
2. Sample container
3. Micro pipette

SPECIMEN COLLECTION

1. Testing should be performed as early as possible after collection. Do not leave serum/Plasma at room temperature for prolonged periods.
2. The samples may be stored at 2°C to 8°C for up to 7 days or frozen at -20°C or lower for up to 30 days.

TEST PROCEDURE

A) For Dengue NS1 Antigen Test:

1. Allow the Test Device and sample reach to room temperature (20°C to 30°C) before opening the foil pouch.
2. Remove the Test Device and plastic Dropper from the pouch and use the test as early as possible.
3. Add two drops (i.e. approximately 90 µl) of serum or plasma sample in well 'S'.
4. Set the timer for 15 minutes.
5. Read the result at 15 minutes. Do not read the result after 20 minutes.

INTERPRETATION OF RESULTS

For Dengue NS1 Antigen Test:

Negative: If pink-red colored line appears at the control side 'C' only.



Positive: If colored lines appear at the control side 'C' and the test side 'T'.



Invalid: The test should be considered invalid if "C" LINE does not appear.

NOTE: The intensity of the pink red color at the test line side (T) will vary depending upon the concentration of dengue virus NS1 antigen in specimen.

Invalid: The test should be considered invalid if "C" LINE does not appear

NOTE: The intensity of the pink red color in the test line region will vary depending upon the infection.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal valid procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

PERFORMANCE CHARACTERISTICS

Internal Evaluation:

Dengue NS1: In an in-house study, total 275 samples were evaluated for sensitivity and specificity. We found the relative sensitivity was 100 % (i. e. 125/125) and the relative specificity was 100 % (i. e. 150/150). The results are summarized in the following table:

Sample	Total Number of samples tested	Rapid Dengue NS1 antigen Test		Sensitivity (%)	Specificity (%)
		Positive	Negative		
Dengue NS1 antigen positive serum Samples	100	100	0	100	-
Dengue NS1 antigen positive plasma Samples	25	25	0	100	-
Negative Human Serum Samples	100	0	100	-	100
Negative Human Plasma Samples	50	0	50	-	100

Cross reactivity was studied using Malaria Pf antigen positive blood samples, Malaria Pv antigen positive blood samples, no cross reactivity was observed.

REFERENCES

1. Pang L, Kitsutani P, Vorndam V, Nakata M, Ayers T, El MJ, Tom T, Reiter P, Rigau-Perez JG, Hayes JM, Mills K, Napier M, Clark GG, Gubler DJ; Hawaii Dengue Outbreak Investigation Team. Dengue fever, Hawaii, 2001-2002. *Emerg Infect Dis.* 2005; 11(5)
2. Alcon-LePoder S, Sivard P, Drouet MT, Talarmin A, Rice C, Flamand M. Secretion of flaviviral non-structural protein NS1: from diagnosis to pathogenesis.
3. Songee L. ranch and Paul N. Levett. Evaluation of four methods for detection of immunoglobulin M antibodies to dengue virus. *Clin. Diagn. Lab. Immunol.* Vol6 (4) p 555-557, 1999.
4. Seth, J. (1991). standardization & quality assurance. In principle and practice of immunoassay, Ed. C.P. Price & D.J. Newman. Macmillan Publishers, pp.154-189.
5. Clinical Evaluation of a rapid immunochromatographic test for the diagnosis of Dengue Virus Infection, Chew Theng Sang, Lim Siew Hoon, Andrea Cuzzubbo, Peter Devine. *Clinical and Diagnostic Laboratory Immunology*, May 1998, Vol. 5, No. 3 p. 407-409.
6. Innis BL, and Nisalak A, et al: An enzyme-linked immunosorbent assay to characterize dengue infections where denude and Japanese encephalitis co-circulate. *Am. J. Trop. Med. Hygiene.* 1989; 40: 418-427. Songee L. ranch and Paul N. Levett. Evaluation of four methods for detection of immunoglobulin M antibodies to dengue virus. *Clin. Diagn. Lab. Immunol.* Vol. 6(4) p 555-557, 1999.



HELICO INTERNATIONAL

Plot # C-30 & 39, KSSIDC Industrial Area,
Doddaballapur-561203 (Karnataka) India.

Contact Details: Tel: +91- 94490-70475

Email: info@helicointernational.com/marketing@helicointernational.com

Website: www.helicointernational.com